

- Run a 1H Experiment First
 - Run with Narrowest possible sweep width
 - run 1H experiment
 - expand around spectrum
 - click on spectrum with left button
 - define left cursor with middle button
 - define right cursor with middle button
 - display will expand around the cursors.
 - click "utilities", then "sw-sfo1"
 - write down "sw" and "O1"
 - click "return"
 - type "rg"
 - write down number
 - DOT NOT CHANGE VALUE
 - hit carriage return on keyboard
- type "dir" and choose "cosy"
- type "acqu" to go to acquisition window
- type "eda"
 - You will see 2 columns of numbers: F2 and F1
 - F2 is the directly detected dimension (normal 1D) (this is t2 in my notes)
 - F1 is the indirectly detected dimension (this is t1 in my notes)
 - You only need to change a few parameters
 - TD in F2 column
 - this is the number of points that will be collected in t2
 - if this number is too large (like what is normally used in 1D H NMR) the data files will be too large for the computer to handle.
 - try 1k -> 2k (k = 1024)
 - TD in F1 column
 - this is the number experiments performed in 2D or equivalently the number of points in t1.
 - if this is too large, not only will the data set be too large to handle, the experiment will take a long time to do.
 - try 512 -> 1k
 - NS (number of scans). Needs to be a multiple of 8.
 - DS (number of dummy scans: scans without acquisition used to create a steady state magnetization for the experiment) try 0, 2, or 4
 - SW F1 and F2. (sweep widths for 1st and 2nd dimensions) enter the value determined in your 1D H NMR experiment (see above)
 - rg (receiver gain) set to value used in 1D H NMR experiment (see above)
 - O1 (center frequency of spectrum) You will need to scroll down a little to find this parameter. Use the value that you used in the 1D H NMR experiment (see above)
 - click on "save"
- Other parameters
 - D1 (recycle delay) set to $(1-5) \cdot T1$ (longitudinal relaxation time)
 - P1 (H 90 pulse) about 10us
- type "expt"
 - This will give you the experiment time
 - If the experiment will take too long, you can either change NS or TD{F1} (in eda)
 - You should expect the experiment to take 3-12 hours.
- Type "zg" to run the experiment.
- Come back when the experiment is finished.

- Set up processing parameters. Type "edp".
 - $SI\{F2\} = TD\{F2\}$
 - $SI\{F1\} = 2 * TD\{F1\}$
 - $wdw\{F2\} = wdw\{F1\} = GM$
 - To set lb and gb for F1 and F2
 - use my defaults of lb = -1.0 and gb = 0.1
 - or find good values for your experiment
 - click "save"
 - type "rser 1"
 - click "process", then "manuel window adjust"
 - will get a dual display
 - Top: damped FID with window function
 - Bottom: Fourier Transform of this FID
 - click "gm"
 - click "+" or "-" next to "lb" and "gb" to adjust the values
 - adjust "lb" and "gb" until FID goes to zero at end and the spectrum looks good to you
 - write down values of "lb" and "gb"
 - click "return"
 - Click "file", then "recall last", then "last 2D data set"
 - type "edp"
 - enter value of "lb" and "gb" into the proper fields
 - click "save"
- type "xfb" to Fourier Transform both dimensions
- 2D spectrum will be displayed
- Phase 2D
 - click "phase"
 - click "row"
 - Will now get cross hairs in 2D spectrum
 - click on a peak with the middle mouse button
 - click with right mouse button to remove cross hairs
 - can change the chosen slice by clicking "+" and "-"
 - when you have the slice that you want, click "mov:1"
 - will move slice to window 1
 - you can change the height of spectrum with "*2" and "/2" buttons below the "mov:1:2:3" buttons
 - choose 2 other rows and move them to windows 2 and 3 with the "mov:2" and the "mov:3" buttons. Choose the to sample entire spectrum (peaks near both edges and in the middle)
 - phase with "PH0" and "PH1" (remember that the peaks in this experiment are not all positive)
 - click "column" and follow the same procedure (steps 2-4 above) to phase the columns.
 - click "return", then "save and return"
the program will ask: "start xfbp?" (perform phasing in both dimensions)
click "OK"
- Plot 2D
 - Define Plot Region
 - click on the box (upper left of xwinmmr window)
 - outline the region of interest
 - hold right mouse button to draw box around peaks
 - click right button to define region
 - Define Intensity of Spectrum

- click "+/-" twice (until gauge on right side of xwinnmr window shows both red and purple region. This shows both positive and negative peaks)
- click on "*8", "/8", "*2", "/2" until the spectrum intensity is as you want it.
- Open "output", then "Define/Show Plot region", then "According to screen limits".
Program will ask:
 - "change levels?" answer "Y"
 - "Please enter number of pos. levels" enter the number of positive levels that you want. try 5
 - "Please enter number of neg. levels" enter the number of negative levels that you want. try 5
 - "Display contours?" answer "N"
- To see what will be plotted
 - type "view"
 - quit view window by clicking on "quit"
- Plot by typing "plot"

II. NOESY

- Run a 1H Experiment First
 - Run with Narrowest possible sweep width
 - run 1H experiment
 - expand around spectrum
 - click on spectrum with left button
 - define left cursor with middle button
 - define right cursor with middle button
 - display will expand around the cursors.
 - click "utilities", then "sw-sfo1"
 - write down "sw" and "O1"
 - click "return"
 - type "rg"
 - write down number
 - DOT NOT CHANGE VALUE
 - hit carriage return on keyboard
- type "dir" and choose "noesy"
- type "acqu" to go to acquisition window
- type "eda"
 - You will see 2 columns of numbers: F2 and F1
 - F2 is the directly detected dimension (normal 1D) (this is t2 in my notes)
 - F1 is the indirectly detected dimension (this is t1 in my notes)
 - You only need to change a few parameters
 - TD in F2 column
 - this is the number of points that will be collected in t2
 - if this number is too large (like what is normally used in 1D H NMR) the data files will be too large for the computer to handle.
 - try 1k -> 2k (k = 1024)
 - TD in F1 column
 - this is the number experiments performed in 2D or equivalently the number of points in t1.
 - if this is too large, not only will the data set be too large to handle, the experiment will take a long time to do.

- try 512 -> 1k
 - NS (number of scans). Needs to be a multiple of 8.
 - DS (number of dummy scans: scans without acquisition used to create a steady state magnetization for the experiment) try 0, 2, or 4
 - SW F1 and F2. (sweep widths for 1st and 2nd dimensions) enter the value determined in your 1D H NMR experiment (see above)
 - rg (receiver gain) set to value used in 1D H NMR experiment (see above)
 - O1 (center frequency of spectrum) You will need to scroll down a little to find this parameter. Use the value that you used in the 1D H NMR experiment (see above)
 - click on "save"
- Other parameters
 - D1 (recycle delay) set to $(1-5)*T1$ (longitudinal relaxation time)
 - P1 (H 90 pulse) about 10us
 - D8 (mixing time for cross relaxation) try $(0.01->1.0)*T1$ (longitudinal relaxation time) you might need to do a series of NOESY experiment as a function of D8
- type "expt"
 - This will give you the experiment time
 - If the experiment will take too long, you can either change NS or TD{F1} (in eda)
 - You should expect the experiment to take 3-12 hours.
- Type "zg" to run the experiment.
- Come back when the experiment is finished.
- Set up processing parameters. type "edp"
 - SI{F2} = TD{F2}
 - SI{F1} = 2*TD{F1}
 - wdw{F2} = wdw{F1} = GM
 - To set lb and gb for F1 and F2
 - use my defaults of lb = -1.0 and gb = 0.1
 - or find good values for your experiment
 - click "save"
 - type "rser 1"
 - click "process", then "manuel window adjust"
 - will get a dual display
 - Top: damped FID with window function
 - Bottom: Fourier Transform of this FID
 - click "gm"
 - click "+" or "-" next to "lb" and "gb" to adjust the values
 - adjust "lb" and "gb" until FID goes to zero at end and the spectrum looks good to you
 - write down values of "lb" and "gb"
 - click "return"
 - Click "file", then "recall last", then "last 2D data set"
 - type "edp"
 - enter value of "lb" and "gb" into the proper fields
 - click "save"
- type "xfb" to Fourier Transform both dimensions
- 2D spectrum will be displayed
- Phase 2D
 - click "phase"
 - click "row"
 - Will now get cross hairs in 2D spectrum
 - click on a peak with the middle mouse button

- click with right mouse button to remove cross hairs
 - can change the chosen slice by clicking "+" and "-"
 - when you have the slice that you want, click "mov:1"
 - will move slice to window 1
 - you can change the height of spectrum with "*2" and "/2" buttons below the "mov:1:2:3" buttons
- choose 2 other rows and move them to windows 2 and 3 with the "mov:2" and "mov:3" buttons. Choose the slices to sample entire spectrum (peaks near both edges and in the middle)
- phase with "PH0" and "PH1" (remember that the peaks in this experiment are not all positive)
- click "column" and follow the same procedure (steps 2-4 above) to phase the columns.
- click "return", then "save and return"
the program will ask: "start xfbp?" (perform phasing in both dimensions)
click "OK"
- Plot 2D
 - Define Plot Region
 - click on the box (upper left of xwinmr window)
 - outline the region of interest
 - hold right mouse button to draw box around peaks
 - click right button to define region
 - Define Intensity of Spectrum
 - click "+/-" twice (until gauge on right side of xwinmr window shows both red and purple region. This shows both positive and negative peaks)
 - click on "*8", "/8", "*2", "/2" until the spectrum intensity is as you want it.
 - Open "output", then "Define/Show Plot region", then "According to screen limits".
Program will ask
 - "change levels?" answer "Y"
 - "Please enter number of pos. levels" enter the number of positive levels that you want. try 5
 - "Please enter number of neg. levels" enter the number of negative levels that you want. try 5
 - "Display contours?" answer "N"
 - To see what will be plotted
 - type "view"
 - quit view window by clicking on "quit"
 - Plot by typing "plot"

III. H C 2D Correlation Experiment (2D HC INEPT)

- Run a 1H and 13C Experiments First
 - Run with both spectra with the narrowest possible sweep widths
 - run 1H experiment
 - expand around spectrum
 - click on spectrum with left button
 - define left cursor with middle button
 - define right cursor with middle button
 - display will expand around the cursors.
 - click "utilities", then "sw-sfo1"
 - write down "sw" and "O1"
 - click "return"

- repeat steps b and c with ^{13}C experiment
- type "dir" and choose "HC_inept_2d"
- type "acqu" to go to acquisition window
- type "eda"
 - You will see 2 columns of numbers: F2 and F1
 - F2 is the directly detected dimension (normal 1D) (this is t2 in my notes)
 - F1 is the indirectly detected dimension (this is t1 in my notes)
 - You only need to change a few parameters
 - TD in F2 column
 - this is the number of points that will be collected in t2
 - if this number is too large (like what is normally used in 1D H NMR) the data files will be too large for the computer to handle.
 - try 1k \rightarrow 2k (k = 1024)
 - TD in F1 column
 - this is the number experiments performed in 2D or equivalently the number of points in t1.
 - if this is too large, not only will the data set be too large to handle, the experiment will take a long time to do.
 - try 512 \rightarrow 1k
 - NS (number of scans). Needs to be a multiple of 16.
 - DS (number of dummy scans: scans without acquisition used to create a steady state magnetization for the experiment) try 0, 2, or 4
 - SW F1 and F2. (sweep widths for 1st and 2nd dimensions)
 - $\text{SW}\{\text{F2}\}$ = sweep width from your ^{13}C NMR experiment (see above)
 - $\text{SW}\{\text{F1}\}$ = sweep width from your ^1H NMR experiment (see above)
 - O1 (center frequency of ^{13}C spectrum) You will need to scroll down a little to find this parameter. Use the value that you used in the ^{13}C experiment (see above)
 - O2 (center frequency of ^1H spectrum) You will need to scroll down a little to find this parameter. Use the value that you used in the ^1H experiment (see above)
 - click on "save"
- Other parameters
 - D1 (recycle delay) set to $(1-5)*\text{T1}$ (longitudinal relaxation time)
 - P1 (H 90 pulse) about 10us
 - type "rga" to set receiver gain.
- type "expt"
 - This will give you the experiment time
 - If the experiment will take too long, you can either change NS or TD{F1} (in eda)
 - You should expect the experiment to take 3-12 hours. THIS TENDS TO BE A LONG EXPERIMENT.
- Type "zg" to run the experiment.
- Come back when the experiment is finished.
- Set up processing parameters. type "edp"
 - $\text{SI}\{\text{F2}\} = \text{TD}\{\text{F2}\}$
 - $\text{SI}\{\text{F1}\} = 2*\text{TD}\{\text{F1}\}$
 - $\text{wdw}\{\text{F2}\} = \text{EM wdw}\{\text{F1}\} = \text{GM}$
 - To set lb and gb for F1 and F2
 - use my defaults of
 - for F2 lb = 5
 - for F1 lb = -1.0 and gb = 0.1
 - or find good values for the lb{F2}

- click "save"
 - type "rser 1"
 - click "process", then "manuel window adjust"
 - will get a dual display
 - Top: damped FID with window function
 - Bottom: Fourier Transform of this FID
 - click "em"
 - click "+" or "-" next to "lb" to adjust the values
 - adjust "lb" until FID goes to zero at end and the spectrum looks good to you
 - write down values of "lb"
 - click "return"
 - Click "file", then "recall last", then "last 2D data set"
 - type "edp"
 - enter value of "lb" into the proper fields
 - click "save"
 - to find values for lb_{F1} and gb_{f1} take a 1H spectrum with $TD = TD_{F1}$ and use the manuel window adjust to find good values. Enter the values in the proper fields in edp.
- type "xfb" to Fourier Transform both dimensions
- 2D spectrum will be displayed
- Phase 2D
 - click "phase"
 - click "row"
 - Will now get cross hairs in 2D spectrum
 - click on a peak with the middle mouse button
 - click with right mouse button to remove cross hairs
 - can change the chosen slice by clicking "+" and "-"
 - when you have the slice that you want, click "mov:1"
 - will move slice to window 1
 - you can change the height of spectrum with "*2" and "/2" buttons below the "mov:1:2:3" buttons
 - choose 2 other rows and move them to windows 2 and 3. Choose the slices to sample entire spectrum (peaks near both edges and in the middle)
 - phase with "PH0" and "PH1" (remember that the peaks in this experiment are ALL POSITIVE)
 - click "column" and follow the same procedure (steps 2-4 above) to phase the columns.
 - click "return", then "save and return"
the program will ask: "start xfbp?" (perform phasing in both dimensions)
click "OK"
- Plot 2D
 - Define Plot Region
 - click on thebox (upper left of xwinmmr window)
 - outline the region of interest
 - hold right mouse button to draw box around peaks
 - click right button to define region
 - Define Intensity of Spectrum
 - click "+/-" twice (until gauge on right side of xwinmmr window shows both red and purple region. This shows both positive and negative peaks)
 - click on "*8", "/8", "*2", "/2" until the spectrum intensity is as you want it.
 - Open "output", then "Define/Show Plot region", then "According to screen limits".
Program will ask

- "change levels?" answer "Y"
- "Please enter number of pos. levels" enter the number of positive levels that you want. try 5
- "Please enter number of neg. levels" enter the number of negative levels that you want. use 0 (there should be no negative peaks in this experiment)
- "Display contours?" answer "N"
- To see what will be plotted
 - type "view"
 - quit view window by clicking on "quit"
- Plot by typing "plot"