

Chemical Behavior of Platinum, Aluminum Oxide and Graphite Crucibles

The following are critical for **PLATINUM**:

- halogens (Cl₂, F₂, Br₂), aqua regia
- Li₂CO₃, prior to emission of CO₂
- SiC from approx. 1000°C
- PbO, FeCl₃
- Be alloys (begin to evaporate just above the melting point)
- HCl with oxidation agents (e.g. chromic acid, manganates, iron(III) salts, molten salts)
- reducing atmospheres
- Pb, Zn, metals such as Pb, Zn, Sn, Ag, Au, Hg, Li, Na, K, Sb, Bi, Ni, Fe, steel, As, Si
- Se above 320°C (immediate cooling and removal of the sample at the end of the measurement recommended to prevent evaporation of selenium)
- metal oxides with reducing substances such as C, organic compounds or H₂
- oxide in an inert gas atmosphere at higher temperatures (reduction)
- sulfur → roughening the surface
- alkali hydroxides, alkali carbonates, alkali sulfates, alkali cyanides and alkali rhodanides at higher temperatures
- KHSO₄ at higher temperatures
- carbon black or free carbon above 1000°C
- SiO₂ under reducing conditions
- SiC and Si₃N₄ at temperatures above 100°C (release of elementary Si)
- HBr, KCN solution at higher temperatures
- high-temperature resistant oxides above 1000°C

no resistance to:

- mixtures of KNO₃ and NaOH at 700°C under exclusion of air
- mixtures of KOH and K₂S at 700°C under exclusion of air
- LiCl at 600°C
- Na₂O₂ at 500°C under exclusion of air
- MgCl₂, Ba(NO₃)₂ at 700°C
- HBr, HJ, H₂O₂ (30%) and HNO₃ at 100°C
- KCl (the decomposition products which form during melting are damaging, melting point: 768°C)

limited resistance to:

- KHF₂, LiF, NaCl at 900°C
- mixtures of NaOH and NaNO₃ at 700°C under exclusion of air

The following are critical for **ALUMINUM OXIDE**:

- N₂ in the presence of carbon: formation of AlN, therefore it is dangerous to measure carbon black at higher temperatures in a nitrogen atmosphere
- F₂: formation of AlF₃ and O₂
- Cl₂: formation of AlCl₃ above 700°C
- sulfur: no reaction with liquid sulfur, in the presence of carbon in the gas phase, formation of sulfides at higher temperatures
- H₂S: formation of up to 3% Al₂S₃ when heated
- C: formation of carbides and Al from approx. 1400°C
- HF: quantitative reaction to Al⁻ and H₂O at higher temperatures
- CuSO₄: diffusion through the bottom of the crucible from approx. 1000°C
- compounds containing silicium, e.g. MoSi₂: contamination of Al₂O₃ from approx. 1200°C in inert gas; reactions at the contact points in air atmosphere
- metal fluorides: attack by the melt, formation of [AlF₆]³⁻ anions and salts similar to cryolites
- SiO₂ glass: glass melts dissolve Al₂O₃
- hydrogen sulfates of alkaline metals and alkaline-earth metals
- HCl: no reaction to 600°C, increased reaction in the presence of carbon at higher temperatures
- B₂O₃ or Borax: melt dissolves Al₂O₃, formation of aluminum borates and aluminum borides
- alkaline and alkaline-earth oxides and their salts with volatile anions: melts form aluminates or double oxides, important, for example, for hydroxides, nitrides, nitrates, carbonates, peroxides, etc.
- CaC₂: formation of Al₄C₃ when heated
- PbO: reaction from 700°C, also important for higher lead oxides and lead salts with anions of volatile acids
- UO₃: reaction begins at 450°C, analogous to PbO
- MeⁿO: Me = Fe²⁺, Cu²⁺, Ni²⁺, etc., formation of spinels; CaO above approx. 1200°C
- alkaline and alkaline-earth ferrites: melt dissolves Al₂O₃
 - zirconium alloys with melting range between 800°C and 1200°C: slow, weak reaction
 - titanium alloys (very high oxygen affinity)
 - reaction with iron-nickel and magnesium alloys possible at higher temperatures

The following are critical for **GRAPHITE**:

- O₂, reaction above 400°C
- molten metals, such as Fe, Co, Ni, Na
- N₂, reaction from 1700°C (formation of small amounts of cyanides)
- oxides (probable reduction upon direct contact)
- water vapor
- F₂, Br₂ at room temperature
- sulfur
- Si at approx. 1400°C (formation of SiC)
- chromic acid (aqueous)
- chlorosulfonic acid ClSO₃H

- SiO₂ → formation of SiC via intermediate product SiO (technical production of SiC above 1800°C; i.e. the reaction between SiO₂ and C definitely begins at lower temperatures)
- nitrous gases (NO, NO₂)
- sulfuric acid, H₂SO₄, concentrated at approx. 150°C, fuming H₂SO₄ at room temperature
- nitric acid, HNO₃, diluted at approx. 90°C, fuming HNO₃ at room temperature
- SO₃, from approx. 100°C
- danger of explosion with perchloric acid, HClO₄
- NaOCl, from approx. 50°C

No claim is made that this overview is exhaustive.

The temperatures given are primarily literature values. Therefore, the temperatures can shift to lower values under experimental conditions. In any case, it is advisable to run tests in separate furnaces.

Compatibility of Calibration Substances with Various Crucible Materials

The table below corresponds to our most current information. No claim is made that it is exhaustive. The user should regard it only as a guide.

Because this is a compilation of data from the literature (H. K. Cammenga et al., Thermochemica Acta, 219, 1993) and from our own experience, it is not possible to cover all possible combinations and measurement conditions. In case of doubt, it is advisable to run preliminary tests in a separate furnace. For platinum crucibles with Al_2O_3 liners, basically the same data like for corundum crucibles applies. The only exception is nickel as a calibration substance. The melting point of nickel with $1455^\circ C$ is above the recommended maximum working temperature for Pt crucibles on a type S sample carrier.

Crucible Material	Calibration Substance																						
	Cyclopentane	Water	Gallium	Indium	Tin	Lead	Zinc	Lithium Sulfate ²⁾	Aluminum	Silver ³⁾	Gold ³⁾	Barium Carbonate ²⁾	Potassium Perchlorate ²⁾	Potassium Chromate ²⁾	Silver Sulfate ²⁾	Rubidium Nitrate ²⁾	Nickel	Bismuth	Cyclohexane	Cesium Chloride ²⁾	Mercury		
Corundum (Al_2O_3)	-	-	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	-	+	-	+	+
Graphite	-	-	+	+	+	+	+	+	+	+	-	?	?	?	?	?	-	?	-	?	?		
Fused Silica	+	+	+	+	+	+	+	-	+	+	-	-	?	?	?	?	-	?	?	?	+		
Aluminum ¹⁾	+	+	-	+	+	+	+	-	-	-	-	-	+	-	+	+	-	+	+	+	+	+	
Silver	+	+	-	-	-	-	-	?	-	-	-	-	+	+	?	+	-	-	?	+	-		
Gold	+	+	-	-	-	-	+	-	-	-	-	?	+	?	?	+	-	-	?	+	-		
Stainless Steel	+	+	-	+	-	+	-	?	-	-	-	-	?	-	-	?	-	+	?	+	+		
Platinum	+	+	-	-	-	-	+	-	-	-	-	+	+	+	+	+	-	-	?	+	-		
Tungsten	-	-	-	?	?	-	+	?	-	+	+	?	?	?	?	?	?	+	?	?	+		

- + No solubility or effect on the melting temperature to be expected.
- ! Corrosion processes possible with negligible changes in the melting temperature.
- - Crucible melts.
 - Melt or transformation product reacts with the crucible material, resulting in greater changes in the melting temperature.
 - Combination not practicable.
- Compatibility unknown.

¹⁾ In air, aluminum forms a thin oxide layer on the surface, which improves the corrosion resistance of the metal. Similar is true for tungsten.

²⁾ Solid-solid transformation.

³⁾ Gold and Silver stick in Al_2O_3 crucibles, but can normally be removed mechanically.