

Sample Preparation. To be able to get absolute MW with dn/dc calculation - concentration of the polymer must be perfectly known. Sample must be clean, dry, and you need to prepare solution with exact concentration (3-4 mg/mL). Prepare solution with THF. As always, make sure your sample is completely dissolved, so you will not filter off part of the sample during filtration. Filter through 0.45 μ m filter, and transfer sample into the autosampler vial (fill between the second and third mark).

- Solvent - THF
- Columns: Two Agilent PLgel MiniMixed-D bed columns + guard (MW linear range: 200-400,000 g/mol)

Use FBS to log into MALS computer & GPC MALS User account on the MALS computer

Loading samples on the Waters 2695 pump Using Empower:

1. STARTING EMPOWER SOFTWARE
 - a. Open Empower icon on desktop (if you need to log in, the log in information is in the folder next to the computer).
 - b. Click on Run Samples
 - c. Open GPC MALS → GPC MALS. Click on Ok.
 - d. Load previously created sample set (MALS mini no UV)
2. Delete unused rows. Insert more rows if necessary- change name, vial #. Load samples into autosampler carousel, and double check vial # in software with position in autosampler carousel.
3. Don't change rows, injection volume or # of injections:
 - i. Rows before samples: Purge Injector; Equilibrate
 - ii. Row after samples: Condition Column
- a. Injection volume: 40 μ L and # of injections: 1.
4. Save, but **do not hit the run (Green Icon) until the Astra software is started!!!!** Note: You may have to enter a date after the Sample Set name (MALS mini no UV) to be able to save the sequence.

Entering samples Using Astra:

1. Open Astra from toolbar (blue icon)
2. Check instrument status by System → Instruments.
 - a. A pop-up window will open. You should see 4 instruments present with green dots.
 - b. You may need to click refresh in the pop-up window. If you do not see 4 instruments with green dots, you may need to restart the computer.
3. Create Astra sequence by File → New → Sequence from Empower
 - a. Log into Empower
 - b. Import Empower Sequence
 - c. Open GPC MALS folder → MALS mini no UV
4. Go to Configuration, Samples and Collection Window:
 - Configuration- Click on the ... → User → GPC MALS mini no UV
 - o Click Apply at bottom of page
5. Samples- Insert sample concentration and dn/dc if known or 0.185 -- for PS (you can change value and recalculate for the correct value of dn/dc later).
 - Do not change run time = 35 min or injection = 40 μ L.
6. Start run with Astra - using green RUN button.
7. Sequence will validate, and you will be able to save your sequence in your folder. After saving, you should see in the Collection window - Waiting for injection.
8. Press start (Green button) in the Empower software
9. Astra will automatically start data collection after sample injection.

RESULTS: Analysis - using ASTRA.

1. Open ASTRA.
2. Open File → Experiment → Sample name.
3. Open Procedures.
 - a. Open Baselines
 - i. Autofind baselines and adjust if necessary → set all
 - ii. Check baseline for each detector by clicking on the checkbox next to the detector number or letter code, adjust if necessary → APPLY
 - iii. When done with baselines → OK.
 - b. Open Peaks- Autofind peaks, can adjust peak window as necessary (To zoom in on baseline, use Ctrl + mouse left click) → APPLY → OK
4. Open Results → Report (summary) to see Mn, Mw, Mp, PDI, etc.
5. Save

Dn/dc value – on line.

NOTE: You must assume 100% injection for calculation of dn/dc.

1. Open Astra.
2. Open File → Experiment → Sample name, unless already open. Right mouse click on experiment → Apply Methods → System → Methods → RI Measurements → 100% Mass Recovery Methods → dn/dc from peak (should open new experiment file on the left side of Astra).
3. Open Procedures.
 - a. Open Baselines, and set baselines if have not done so already
 - i. Autofind baselines and adjust if necessary → set all
 - ii. Check baseline for each detector by clicking on the checkbox next to the detector number or letter code, adjust if necessary → APPLY
 - iii. When done with baselines → OK.
 - b. Open Peaks
 - i. Autofind peaks, can adjust peak window as necessary (To zoom in on baseline, use Ctrl + mouse left click)
 - ii. Check concentration of injected mass → APPLY → OK
4. Open Results → Report (summary) New dn/dc value is under peak results
5. Enter new dn/dc value in Source Experiment file → Procedures → Peaks, type in dn/dc value
6. Open Results of Source Experiment File → Report (summary) to see recalculated Mn, Mw, MP
7. For completely new polymers or difficult systems it is possible to calculate dn/dc for 5 perfect concentrations of the polymer (for same solvent, temperature and wavelength of the laser)

Shutting down:

1. Check that system changes from waste to recycle mode.
2. Turn off laser on Heleos if not done by software